


The Golden Gate™ - Low temp. Diamond ATR System

A variable temperature diamond ATR system for the FTIR measurement of solids and liquids at ambient to near liquid nitrogen temperatures.



Key Sampling Benefits

- Continuous control of sample temperature from near liquid nitrogen temperatures to 80 °C
- High thermal conductivity provides rapid cooling and temperature stabilization
- Proven, strong clamping device, based on the Golden Gate™ Diamond ATR, allows rapid, reproducible sample throughput
- Thermally insulated copper and stainless steel dewar allows for the use of liquid nitrogen, dry ice, or salt and water mixtures

The Golden Gate low temperature diamond ATR system is the first ATR accessory to provide high performance ATR measurements down to near liquid nitrogen temperature.

The system uses a thermally insulated copper and stainless steel dewar in conjunction with an integral heater to provide continuous temperature control across the working range.

The high thermal conductivity of the diamond crystal provides rapid temperature stabilization, accurate temperature measurement, and avoids temperature gradients across the sample.

The diamond ATR crystal is high temperature metal bonded into a tungsten carbide support disk and the top plate is hard-anodized to make the accessory chemically resistant and capable of withstanding the pressures required for

optimum crystal and sample contact. Pressure is applied to the sample using the quick lock and release bridge from the proven Golden Gate diamond ATR system. This clamping device allows rapid sample throughput and reproducible solid sampling. It also has a built-in torque limiter to control the loads applied to the diamond.

The Golden Gate low temperature diamond ATR is constructed in such a way that the crystal mounting is under a constant load. This ensures that the diamond is kept in constant optical alignment, negating the effects of thermal expansion and contraction. The upper body of the system is separated from the top-plate by a replaceable, thermally conducting spacer, avoiding the need for use of sealants that may contaminate the diamond surface. The system is easily and quickly taken apart for cleaning.

The Golden Gate™ - Low Temperature Diamond ATR System




Figure 1: PTFE recorded at room temperature (black) and at -100 °C (blue) using a DTGS detector at 4 cm⁻¹ resolution.




Figure 2: Carbon tetrachloride measured at room temperature (blue) and -143 °C (black) using a DTGS detector at 4 cm⁻¹ resolution.

Ordering information

GS10592

Golden Gate Low Temperature

Diamond ATR System, includes optics unit with choice of ZnSe or KRS -5 lenses, low temperature diamond top plate and temperature controller, sample clamping system with integrated torque limiter, Benchmark™ Instrument Base-plate and purge bellows.

Please specify instrument make and model, required voltage and country of destination.

The temperature controller is supplied without any computer connection, if required, please specify RS232, RS485 or USB connection.

GS10590

Golden Gate™ Low Temperature Diamond ATR Top Plate with temperature controller only, includes sample clamping system with integrated torque limiter.

Please specify instrument make and model, required voltage and country of destination.

Spares and Consumables

GS10594 High thermal transfer gaskets (pkt 5)

Options

GS28000 RS232 Connection kit
 GS28001 USB Connection kit
 GS28002 RS485 Connection kit

• Illustrations, descriptions and specifications in this data sheet were correct at the time of going to press. However, Specac's policy is one of continuous product development and we reserve the right to change descriptions and specifications at any time. For the latest details please contact your local Specac office or representative.

SPECAC LTD.,
 River House, 97 Cray Ave, St. Mary Cray, Orpington,
 Kent, BR5 4HE UK
 T : +44 (0)1689 873134 F : +44 (0)1689 878527
 E : sales@specac.co.uk Registered in England No. 1008689

SPECAC INC.,
 50 Sharpe Drive, Cranston, RI 02920 USA
 T : Toll Free +1 800 447 2558
 E : sales@specac.co.uk W : www.specac.com
 Specac is part of Smiths Group plc

7 - 05/11