

Press Family

Atlas™ Heated Platens

A quick and simple accessory for converting standard Hydraulic Presses to heatable versions.

The Atlas™ Heated Platens have been designed to provide heatable pressing surfaces compatible with all Specac Hydraulic Presses. The platens are easily installed by replacing the bolster fixing of the press by the top platen and placing the bottom platen over the lower pressing face of the press.

The platens have a large surface area, are rugged, durable and have full CE safety approval. They are fitted with a permanent thermocouple which monitors the temperature close to the pressing surface. The press connection ports are water cooled which allows for an efficient heating and cooling cycle.

The platens are controlled by a digital automatic temperature controller with a stability of $\pm 1^\circ\text{C}$.

Key Features

- Digital temperature display
- Controlled temperature to 300°C with water cooling
- 15 Ton load bearing capacity
- CE Safety Approved
- 100mm diameter pressing surface
- Compatible with all Specac Atlas™ Hydraulic Presses



Ordering Information

GS15515 Atlas™ Water Colled Heated Platens and Automatic Temperature Controller with Digital Display.
Includes 2 meters of 6mm id. PVC tubing and flow connectors.

Specify 220V or 110V and country of usage

Spares and Consumables

GS15601 PTFE sheets 0.25mm thick for Water Cooled Heated Platens (10)
GS15512 Replacement heater set for Platens (220V)
GS15522 Replacement heater set for Platens (110V)

Illustrations, descriptions and specifications in this datasheet were correct at the time of going to press. However, Specac's policy is one of continuous product development and we reserve the right to change descriptions and specifications at any time.

For the latest details please contact your local Specac office or representative.

SPECAC LTD.,
River House 97 Cray Avenue, Orpington,
Kent. BR5 4HE UK
T: +44 (0) 1689 873134 F: +44 (0) 1689 878527
E: sales@specac.co.uk Registered in England No. 1008689

SPECAC INC.
414 Commerce Drive, #175
Fort Washington, PA 19034, USA
T: +1 866 726 1126
E: sales@specac.com