

Window Polishing Kit

Keep windows in top condition with the Specac Window Polishing Kit



The window polishing kit contains all the essential materials required to clean and repolish infrared transmission windows.

Repolishing can be achieved efficiently and economically with a minimum degree of skill.

Full instructions are included and all consumable parts are replaceable.

• Illustrations, descriptions and specifications in this data sheet were correct at the time of going to press. However, Specac's policy is one of continuous product development and we reserve the right to change descriptions and specifications at any time. For the latest details please contact your local Specac office or representative.

Ordering Information

GS04000 Window Polishing Kit

The Window Polishing Kit includes the following items. (Each item can also be ordered separately as spares)

GS04010 (5) spare Polishing Selvyts (self adhesive pads for Glass polishing lap)

GS04090 (5) spare lapping papers (for smoothing lap)

GS04040 Polyethylene bottle (2)

GS04050 Cleaning sponge

GS04060 Bottle of polishing rouge (approx. 75g net)

GS04080 Non adhesive plain Selvyt strips (2)

GS04095 Brushes (2)

SPECAC LTD.,
River House, 97 Cray Ave, St. Mary Cray, Orpington,
Kent, BR5 4HE UK
T : +44 (0)1689 873134 F : +44 (0)1689 878527
E : sales@specac.co.uk Registered in England No. 1008689

SPECAC INC.,
414 Commerce Drive, #175
Fort Washington, PA
19034, U.S.A.
T : 1 866 726 1126
E : sales@specac.com