

# SpecacDie™



**For producing KBr pellets without using a hydraulic press**

## Applications

- Halide pellets for IR transmission analysis

## Features

- Integrated pellet holder
- Rapid sample preparation
- Evacuatable for sample clarity and quality
- 8.5mm diameter pellets

## Ordering Information

- 03800** SpecacDie™ with bolts and seals only
- 03700** SpecacDie™ Kit  
Includes: SpecacDie™, SpecacDie™ Holder, bench mounted wrench, open ended wrench, spare set of bolts, spare set of seals, bottle of KBr powder (50 g)
- 03560** Circular SpecacDie™ holder
- 03570** Set of bolts
- 03590** Set of two open ended wrenches, 9/16" (14.3 mm) and 1/2" (12.7 mm) AF
- 03595** Bench mounted wrench
- 03580** Set of seals
- 03610** KBr powder

The SpecacDie™ is an easy to use halide pellet maker ideal for the laboratory with small numbers of samples to analyze.

Powdered samples are loaded into the SpecacDie™ body and pressed into a thin disk between two highly polished bolt ends.

Once the bolts are removed, the die body acts as the pellet holder which is then mounted on the 3" x 2" holder within the spectrometer for analysis.

THE SOLUTION IS SPECAC WHAT'S THE PROBLEM?

smiths  
A part of Smiths Group plc

# Specac

www.specac.com



**SPECAC LIMITED,**  
River House, 97 Cray Ave, St. Mary Cray,  
Orpington, Kent BR5 4HE U.K.  
T : +44 (0)1689 873134  
F : +44 (0)1689 878527  
E : sales@specac.com  
Registered No. 1008689 England

**SPECAC INC.,**  
50 Sharpe Drive,  
Cranston, RI 02920 U.S.A.  
T : Toll Free +1 800 447 2558  
E : sales@specac.com  
W : www.specac.com  
Specac is part of Smiths Group plc