

### The Minidiff Plus Diffuse Reflectance Accessory For FT-IR Spectroscopy



The Minidiff Plus diffuse reflectance accessory is designed for routine sampling. High performance solids analysis is made simple through the optimised zero alignment optics and Benchmark™ baseplate mounting system. Sampling versatility is maintained through the choice of either standard solids cup or Diabrase™ abrasive sampling pads. Up to 3 samples can be loaded at the same time, while the sample introduction system reduces the risk of sample spillage.

## Features

### Zero Alignment

- Completely pre-aligned with simple sample introduction/removal.

### Optimised Throughout

- Designed for a high throughput on all FT-IR spectrometers.

### Maximum Versatility

- Sample cups for powders and 'Diabrase' sampler for non particulate solid samples.

### Budget Cost

- Cost engineered to meet the budget of all laboratories.

## Ordering Information

**4510** Minidiff Plus diffuse reflectance accessory.  
Includes: Optics system with Benchmark™ instrument baseplate, abrasive sample holders (2), sample cup holders (2), Diabrase™ abrasive pads (20), KBr powder (50g) and pestle and mortar.

*Please specify instrument make and model.*

**4505** Spare sample cup holders with 3 cups (2)  
**4506** Spare packet of Diabrase™ sampling pads (100)  
**4508** Spare abrasive sample holders with 3 pads (2)  
**3610** KBr powder (50g)  
**3600** Pestle and mortar

THE SOLUTION IS SPECAC WHAT'S THE PROBLEM?

smiths  
A part of Smiths Group plc

# Specac

www.specac.com



**SPECAC LIMITED,**  
River House, 97 Cray Ave, St. Mary Cray,  
Orpington, Kent BR5 4HE U.K.  
T : +44 (0)1689 873134  
F : +44 (0)1689 878527  
E : sales@specac.com  
Registered No. 1008689 England

**SPECAC INC.,**  
50 Sharpe Drive,  
Cranston, RI 02920 U.S.A.  
T : Toll Free +1 800 447 2558  
E : sales@specac.com  
W : www.specac.com  
Specac is part of Smiths Group plc