

Sample Storage and Holders




Specacard™ for mounting thin polymer films and KBr pellets for analysis and storage

.....

Features

- Self sealing adhesive coating
- Indexing space
- Clear aperture: 10 mm diameter or 25x10 mm
- Slide mounted (3"x2")
- Filing rack available



Magnetic Film Holder for mounting thin polymer films for analysis

.....

Features

- Magnetic retarding ring
- Clear aperture: 10 mm diameter
- Slide mounted (3"x2")

.....

Ordering Information

GS03800 Specacards 10 mm dia. aperture (100)

GS03810 Specacards 25x10 mm aperture (100)

GS03815 Filing rack for Specacards

GS03820 Magnetic Film Holder

• Illustrations, descriptions and specifications in this data sheet were correct at the time of going to press. However, Specac's policy is one of continuous product development and we reserve the right to change descriptions and specifications at any time. For the latest details please contact your local Specac office or representative.

Sample Storage and Holders



Specacabinet™ for storage of hygroscopic IR materials under controlled temperature

Features

- Thermostatically controlled from ambient to 80°
- Size (L x W x H): Internal: 737 x 331 x 432mm
External: 762 x 356 x 508mm
- Can be wall mounted

Ordering Information

- GS19100** Specacabinet (220 V, 50Hz) UK
- GS19102** Specacabinet (220 V, 50 Hz) EU

SPECAC LTD.,
River House, 97 Cray Ave, St. Mary Cray, Orpington,
Kent, BR5 4HE UK
T : +44 (0)1689 873134 F : +44 (0)1689 878527
E : sales@specac.co.uk Registered in England No. 1008689

SPECAC INC.,
50 Sharpe Drive, Cranston, RI 02920 USA
T : Toll Free +1 800 447 2558
E : sales@specac.co.uk W : www.specac.com
Specac is part of Smiths Group plc