

The Minidiff Plus - Diffuse Reflectance Accessory For FT-IR Spectroscopy

The Minidiff Plus diffuse reflectance accessory is designed for routine sampling. High performance solids analysis is made simple through the optimised zero alignment optics and Benchmark™ baseplate mounting system. Sampling versatility is maintained through the choice of either standard solids cup or Diabrase™ abrasive sampling pads. Up to 3 samples can be loaded at the same time, while the sample introduction system reduces the risk of sample spillage.



• Illustrations, descriptions and specifications in this data sheet were correct at the time of going to press. However, Specac's policy is one of continuous product development and we reserve the right to change descriptions and specifications at any time. For the latest details please contact your local Specac office or representative.

.....
Features

Zero Alignment

- Completely pre-aligned with simple sample introduction/removal.

Optimised Throughout

- Designed for a high throughput on all FT-IR spectrometers.

Maximum Versatility

- Sample cups for powders and 'Diabrase' sampler for non particulate solid samples.

Budget Cost

- Cost engineered to meet the budget of all laboratories.

.....
Ordering Information

GS4510 Minidiff Plus diffuse reflectance accessory.

Includes: Optics system with Benchmark™ instrument baseplate, abrasive sample holders (2), sample cup holders (2), Diabrase™ abrasive pads (20), KBr powder (50g) and pestle and mortar.

Please specify instrument make and model.

GS4505 Spare sample cup holders with 3 cups (2)

GS4506 Spare packet of Diabrase™ sampling pads (100)

GS4508 Spare abrasive sample holders with 3 pads (2)

GS3610 KBr powder (50g)

GS3600 Pestle and mortar

SPECAC LTD.,
River House, 97 Cray Ave, St. Mary Cray, Orpington,
Kent, BR5 4HE UK
T : +44 (0)1689 873134 F : +44 (0)1689 878527
E : sales@specac.co.uk Registered in England No. 1008689

SPECAC INC.,
50 Sharpe Drive, Cranston, RI 02920 USA
T : Toll Free +1 800 447 2558
E : sales@specac.co.uk W : www.specac.com
Specac is part of Smiths Group plc

9 - 05/11