

Sample Prep**Low Tonnage Gauge Kit**

Low Tonnage Conversion Kit for Atlas 15T Manual Hydraulic Press

The low tonnage conversion kit enables a standard 15T Atlas Manual Hydraulic Press to be used more accurately for reading of a low tonnage load that has been applied during a pressing, due to its finer divisions. The low tonnage gauge conversion kits provide an additional low tonnage load gauge to be used along with the standard 15 ton load gauge fitted on the press.

The appropriate low tonnage gauge kit is fitted by the user to allow both gauges to be connected to the press and used in-line to register an applied tonnage load. The lower tonnage gauge that has been fitted in this way to a press can be independently switched off from the pressure system, but any pressure of oil in the system will always be registered as an appropriate applied tonnage load at the 15 ton gauge fitted.

Key Features

- Easy to Install and Use
- Accurate Low Tonnage Display



Ordering Information

GS15051	Gauge Conversion Kit 0-1 Ton
GS15052	Gauge Conversion Kit 0-2 Ton
GS15055	Gauge Conversion Kit 0-5 Ton

Illustrations, descriptions and specifications in this datasheet were correct at the time of going to press. However, Specac's policy is one of continuous product development and we reserve the right to change descriptions and specifications at any time.

For the latest details please contact your local Specac office or representative.

SPECAC LTD.,
River House 97 Cray Avenue, Orpington,
Kent. BR5 4HE UK
T: +44 (0) 1689 873134 F: +44 (0) 1689 878527
E: sales@specac.co.uk Registered in England No. 1008689

SPECAC INC.
414 Commerce Drive, #175
Fort Washington, PA 19034, USA
T: +1 866 726 1126
E: sales@specac.com